



CooperVision®  
Live Brightly.®

# Product Specifications Guide



## MyDay® daily disposable family

### MyDay® daily disposable



High-performing 1 Day contact lenses for demanding lifestyles.

#### Wearing Schedule

Daily disposable

#### Carton Size Options

30 or 90 pack

#### Material / H<sub>2</sub>O Content

stenfilcon A / 54%

**Base Curve** 8.4mm

**Diameter** 14.2mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** Yes

#### Sphere Power

+8.00 to -12.00  
(to +5.00 & -6.00D  
in -0.25D Steps)  
(from +5.00 &  
-6.00 in 0.50D  
Steps)

#### Design

Aspheric

### MyDay® daily disposable toric



High-performing 1 Day contact lenses for wearers with astigmatism.

#### Wearing Schedule

Daily disposable

#### Carton Size Options

30 pack

#### Material / H<sub>2</sub>O Content

stenfilcon A / 54%

**Base Curve** 8.6mm

**Diameter** 14.5mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** Yes

#### Sphere Power

-0.00D to -6.00D  
(0.25D Steps)

#### Cylinder Power

-0.75  
-1.25  
-1.75

#### Axis

10°, 20°, 30°,  
40°, 50°, 60°,  
70°, 80°, 90°,  
100°, 110°, 120°,  
130°, 140°, 150°  
160°, 170°, 180°

#### Design

Optimized  
Toric Lens  
Geometry™

-2.25

10°, 20°, 90°,  
160°, 170°, 180°

-6.50 to -10.00  
(0.50D Steps)

-0.75  
-1.25  
-1.75

10°, 20°, 70°,  
80°, 90°, 100°,  
110°, 160°, 170°,  
180°

+ 0.25D to +6.00D  
(0.25D Steps)

-6.50 to -10.00  
(0.50D Steps)

-2.25

10°, 20°, 90°,  
160°, 170°, 180°

+0.50D to +6.00D  
(0.50D Steps)

## clariti<sup>®</sup> 1 day family

### clariti<sup>®</sup> 1 day



Great value and high breathability<sup>1</sup>  
in a convenient daily disposable  
contact lens.

#### Wearing Schedule

Daily disposable

#### Carton Size Options

30 or 90 pack

#### Material / H<sub>2</sub>O Content

somofilcon A / 56%

**Base Curve** 8.6mm

**Diameter** 14.1mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** No

#### Sphere Power

+8.00 to -10.00  
(to +/- 6.00 in 0.25  
steps)  
(from +/- 6.00 in  
0.50 steps)

#### Design

Aspheric

### clariti<sup>®</sup> 1 day toric



Great value and high breathability<sup>1</sup>  
in a convenient daily disposable  
contact lens.

#### Wearing Schedule

Daily disposable

#### Carton Size Options

30 pack

#### Material / H<sub>2</sub>O Content

somofilcon A / 56%

**Base Curve** 8.6mm

**Diameter** 14.3mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** No

#### Sphere Power

Plano to -9.00D  
(above -6.00 in  
0.50 steps)

#### Cylinder Power

-0.75  
-1.25  
-1.75

#### Axis

10°, 20°, 60°,  
70°, 80°, 90°,  
100°, 110°, 120°,  
160°, 170°, 180°

Plano to -9.00D  
(above -6.00 in  
0.50 steps)

-2.25

10°, 20°, 90°,  
160°, 170°, 180°

+0.25 to +4.00  
(0.25 steps)

-0.75  
-1.25  
-1.75

20°, 70°, 90°,  
110°, 160°, 180°

**Design**

Aspheric

## clariti® 1 day multifocal



Great value and high breathability<sup>1</sup> in the convenience of a daily disposable contact lens.

**Wearing Schedule**  
Daily disposable

**Carton Size Options**  
30 pack

**Material / H<sub>2</sub>O Content**  
somofilcon A / 56%

**Base Curve** 8.6mm

**Diameter** 14.1mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** No

**Sphere Power**  
+5.00D to -6.00D  
(in 0.25 Steps)

**Add Power**  
Low:  
up to +2.25D

High:  
+2.50D to  
+3.00D

**Design**  
Centre near,  
progressive  
intermediate, and  
peripheral distance

## Biofinity® family

### Biofinity Energys®



Technologically advanced contact lenses for your digital lifestyle.

**Wearing Schedule**  
Daily or up to 29  
days continuous

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
comfilcon A / 48%

**Base Curve** 8.6mm

**Diameter** 14.0mm

**FDA Group** SiHy

**UV Blocking\*** No

**Visibility Tint** Yes

**Sphere Power**  
+8.00D to -12.00D  
(0.50D steps after  
+/-6.00D)

**Design**  
Aspheric  
Digital Zone  
Optics® lens  
design

### Biofinity®



Extraordinary visual performance.  
Premium all-day comfort.

**Wearing Schedule**  
Daily or up to 29  
days continuous

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
comfilcon A / 48%

**Base Curve** 8.6mm

**Diameter** 14.0mm

**FDA Group** SiHy

**UV Blocking\*** No


**Visibility Tint** Yes

**Sphere Power**  
+8.00D to -12.00D  
(0.50D steps after  
+/-6.00D)

**Design**  
Aspheric

+15.00D to  
-20.00D  
(0.50D steps after  
+/-6.00D)

## Biofinity® toric



Extraordinary visual performance.  
Premium all-day comfort.

**Wearing Schedule**  
Daily or up to 29  
days continuous

**Carton Size  
Options**  
3 pack

**Material / H<sub>2</sub>O  
Content**  
comfilcon A / 48%

**Base Curve** 8.7mm  
**Diameter** 14.5mm  
**FDA Group** SiHy

**Sphere Power**  
+8.00D to -10.00D  
(0.50D steps after  
+/-6.00D)

**Cylinder  
Power**  
-0.75  
-1.25  
-1.75  
-2.25

**Axis**  
10° to 180°  
(10° steps)

**UV Blocking\*** No  
**Visibility Tint** Yes

**Design**  
Optimized  
Toric Lens  
Geometry™

+8.50 to +10.00  
(0.50D steps)

-0.75 to -2.25  
(in 0.50 steps)

5° to 180°  
(5° Steps)

+10.00 to -10.00  
(0.50D steps after  
+/- 6.00D)

-2.75 to -5.75  
(in 0.50 steps)

5° to 180°  
(5° Steps)

## Biofinity® multifocal



Extraordinary visual performance  
at every distance. Premium all-day  
comfort.

**Wearing Schedule**  
Daily wear or  
continuous wear up  
to 29 days

**Carton Size  
Options**  
6 pack

**Material / H<sub>2</sub>O  
Content**  
comfilcon A / 48%

**Base Curve** 8.6mm  
**Diameter** 14.0mm  
**FDA Group** SiHy

**Sphere Power**  
+6.00D to -10.00D  
(0.50D steps after  
-6.00D)

**Add Power**  
+1.00  
+1.50  
+2.00  
+2.50

**Design**  
D lens  
N lens

**UV Blocking\*** No  
**Visibility Tint** Yes

## Avaira Vitality™ family

### Avaira Vitality™



High performance. Outstanding  
value.

**Wearing Schedule**  
Monthly

**Carton Size  
Options**  
6 pack

**Material / H<sub>2</sub>O  
Content**  
fanfilcon A / 55%

**Base Curve** 8.4mm  
**Diameter** 14.2mm  
**FDA Group** SiHy

**Sphere Power**  
+08.00 to -12.00  
(to +/- 6.00 in 0.25  
steps)  
(from +/- 6.00 in  
0.50 steps)

**Design**  
Aspheric front  
surface

**UV Blocking\*** Yes  
**Visibility Tint** Yes

## Avaira Vitality™ toric



**Wearing Schedule**  
Monthly

**Carton Size Options**  
3 pack

**Material / H<sub>2</sub>O Content**  
fanfilcon A / 55%

**Base Curve** 8.5mm

**Diameter** 14.5mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** Yes

**Sphere Power**  
-10.00D to +8.00D  
(0.50D steps after  
+/-6.00D)

**Cylinder Power**  
-0.75  
-1.25  
-1.75  
-2.25

**Axis**  
10° to 180°  
(10° steps)

**Design**  
Optimized Toric  
Lens Geometry™

High performance. Outstanding value.

## Proclear® family

### Proclear® 1 day



**Wearing Schedule**  
Daily

**Carton Size Options**  
30 or 90 pack

**Material / H<sub>2</sub>O Content**  
omafilcon A / 60%

**Base Curve** 8.7mm

**Diameter** 14.2mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

**Sphere Power**  
+8.00D to -12.00D  
(0.50D steps  
after +5.00D and  
-6.00D)

**Design**  
Aspheric

Clear vision\*, convenience and comfort, for contact lens wearers who experience eye dryness.

### Proclear® 1 day multifocal



**Wearing Schedule**  
Daily disposable

**Carton Size Options**  
30 pack

**Material / H<sub>2</sub>O Content**  
omafilcon A / 60%

**Base Curve** 8.7mm

**Diameter** 14.2mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

**Sphere Power**  
+6.00D to -10.00D  
(0.50D steps after  
-6.00D)

**Add Power**  
Single power  
profile, up to  
+2.50

**Design**  
Centre near  
aspheric

Convenience and comfort for multifocal contact lens wearers who experience eye dryness.

## Proclear®



Clear vision<sup>2</sup> and comfort for contact lens wearers who experience eye dryness.

<b>Wearing Schedule</b> Daily	<b>Carton Size Options</b> 6 pack	<b>Material / H<sub>2</sub>O Content</b> omafilcon B / 62%	<b>Base Curve</b> 8.6mm
<b>Sphere Power</b> +20.00D to -20.00D (0.50D steps after +/-6.00D)	<b>Design</b> Sphere		<b>Diameter</b> 14.2mm
			<b>FDA Group</b> 2
			<b>UV Blocking*</b> No
			<b>Visibility Tint</b> Yes

## Proclear® toric



Helps provide all-day comfort for contact lens wearers with astigmatism who experience eye dryness.

<b>Wearing Schedule</b> Daily	<b>Carton Size Options</b> 3 or 6 pack	<b>Material / H<sub>2</sub>O Content</b> omafilcon A / 62%	<b>Base Curve</b> 8.8mm
<b>Sphere Power</b> +6.00D to -8.00D (0.50D steps after -6.50D)	<b>Cylinder Power</b> -0.75 -1.25 -1.75 -2.25	<b>Axis</b> 10° to 180° (10° steps)	<b>Diameter</b> 14.4mm
+6.25D to +10.00D (0.50D steps after +6.50D)	-0.75 to -2.25 (0.50D steps)	5° to 180° (5° steps)	<b>FDA Group</b> 2
-8.50D to -10.00D (0.50D steps)			<b>UV Blocking*</b> No
			<b>Visibility Tint</b> Yes
+10.00D to -10.00D (0.50D steps, +6.00D to -6.50D in 0.25D steps)	-2.75 to -5.75 (0.50D steps)	5° to 180° (5° steps)	

## Proclear<sup>®</sup> multifocal



Helps provide all-day comfort for multifocal contact lens wearers who experience eye dryness.

**Wearing Schedule**  
Daily

**Carton Size Options**  
3 or 6 pack

**Material / H<sub>2</sub>O Content**  
omafilcon A / 62%

**Base Curve** 8.7mm

**Diameter** 14.4mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

**Sphere Power**

+6.00 to -08.00

**Add Power**

+1.00  
+1.50  
+2.00  
+2.50

**Design**

D lens  
N lens

+20.00D to  
-20.00D

+3.00  
+3.50  
+4.00

+6.00D to -6.00D  
(0.25D steps)

+6.00D to +20.00D  
(0.50D steps)

-6.50D to -20.00D  
(0.50D steps)

## Proclear<sup>®</sup> multifocal toric



Helps provide all-day comfort for multifocal contact lens wearers with astigmatism who experience eye dryness.

**Wearing Schedule**  
Daily

**Carton Size Options**  
3 pack

**Material / H<sub>2</sub>O Content**  
omafilcon B / 62%

**Base Curve** 8.4, 8.8mm

**Diameter** 14.4mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

**Sphere Power**

+20.00D to  
-20.00D  
(0.50D steps after  
+/-6.50D)

**Cylinder Power**

-0.75 to -5.75  
(0.50D steps)

**Axis**

5° to 180°  
(5° steps)

**Design**

D lens  
N lens

\* Due to aspheric optics.

<sup>1</sup> Brennan NA. Beyond Flux: Total Corneal Oxygen Consumption as an Index of Corneal Oxygen During Contact Lens Wear. Optom Vis Sci 2005. 2005;82(6):467-472.

<sup>2</sup> CVI data on file 2017 & 2014; non-dispensing short-term fitting study. N=40 & n=20.

\* UV-blocking contact lenses help provide protection against transmission of harmful UV radiation to eye but are not substitutes for protective UV-absorbing eyewear, such as UV absorbing goggles or sunglasses, as they do not completely cover the eye or surrounding area. Continue to use UV-absorbing eyewear as directed by your eye care professional.

All trademarks, service marks, trade names, trade dress and product names featured herein are registered trademarks and trademarks of CooperVision, its affiliates and subsidiaries, unless otherwise noted. CooperVision's trademarks may be used publicly only with permission.