




CooperVision®
Live Brightly.™

Product Specifications Guide



MyDay® daily disposable family

MyDay® daily disposable



High-performing 1 Day contact lenses for demanding lifestyles.

Wearing Schedule
Daily disposable

Sphere Power
+8.00 to -12.00
(to +5.00 & -6.00D
in -0.25D Steps)
(from +5.00 &
-6.00 in 0.50D
Steps)

Carton Size
Options
30 or 90 pack

Design
Aspheric

Material / H₂O
Content
stenfilcon A / 54%

Base Curve 8.4mm
Diameter 14.2mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint Yes

MyDay® daily disposable toric



High-performing 1 Day contact lenses for wearers with astigmatism.

Wearing Schedule
Daily disposable

Sphere Power
-0.00 to -6.00
(0.25D Steps)
-6.50 to -10.00
(0.50D Steps)

Carton Size
Options
30 pack
Cylinder
Power

-0.75
-1.25
-1.75
-2.25

Material / H₂O
Content
stenfilcon A / 54%
Axis

Around the clock
10° - 180°

Base Curve 8.6mm
Diameter 14.5mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint Yes
Design
Optimized
Toric Lens
Geometry™

+0.25 to +6.00
(0.25D Steps)
+6.50 to +8.00
(0.50D Steps)

-0.75
-1.25
-1.75
-2.25

Around the clock
10° - 180°

Brilliant Futures™ Myopia Management Program with MiSight® 1 day

Brilliant Futures™ Myopia Management Program with MiSight®



Wearing Schedule

Daily

Sphere Power

-0.50D to -6.00D
(0.25D steps)

Carton Size
Options

90 pack

Material / H₂O
Content

omafilcon A / 60%

Base Curve 8.7mm

Diameter 14.2mm

FDA Group 2

UV Blocking* No

Visibility Tint Yes

Child-friendly daily disposable contact lenses proven to slow the progression of nearsightedness in children aged 8-12 at the initiation of treatment.²⁺

clariti® 1 day family

clariti® 1 day



Great value and high breathability¹ in a convenient daily disposable contact lens.

Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.6mm
Daily disposable	30 or 90 pack	somofilcon A / 56%	Diameter	14.1mm
Sphere Power	Design		FDA Group	SiHy
+8.00 to -10.00 (to +/- 6.00 in 0.25 steps) (from +/- 6.00 in 0.50 steps)	Aspheric		UV Blocking*	Yes
			Visibility Tint	No

clariti® 1 day toric



Great value and high breathability¹ in a convenient daily disposable contact lens.

Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.6mm
Daily disposable	30 pack	somofilcon A / 56%	Diameter	14.3mm
Sphere Power	Cylinder Power	Axis	FDA Group	SiHy
Plano to -9.00D (above -6.00 in 0.50 steps)	-0.75 -1.25 -1.75	10°, 20°, 60°, 70°, 80°, 90°, 100°, 110°, 120°, 160°, 170°, 180°	UV Blocking*	Yes
			Visibility Tint	No
			Design	Aspheric
Plano to -9.00D (above -6.00 in 0.50 steps)	-2.25	10°, 20°, 90°, 160°, 170°, 180°		
+0.25 to +4.00 (0.25 steps)	-0.75 -1.25 -1.75	20°, 70°, 90°, 110°, 160°, 180°		

clariti® 1 day multifocal



Great value and high breathability¹ in the convenience of a daily disposable contact lens.

Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.6mm
Daily disposable	30 pack	somofilcon A / 56%	Diameter	14.1mm
Sphere Power	Add Power	Design	FDA Group	SiHy
+5.00D to -6.00D (in 0.25 Steps)	Low: up to +2.25D High: +2.50D to +3.00D	Centre near, progressive intermediate, and peripheral distance	UV Blocking*	Yes
			Visibility Tint	No

Biofinity® family

Biofinity Energys®



Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.6mm
Daily or up to 29 days continuous	6 pack	comfilcon A / 48%	Diameter	14.0mm
Sphere Power	Design		FDA Group	SiHy
+8.00D to -12.00D (0.50D steps after +/-6.00D)	Aspheric Digital Zone Optics® lens design		UV Blocking*	No
			Visibility Tint	Yes

Technologically advanced contact lenses for your digital lifestyle.

Biofinity®



Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.6mm
Daily or up to 29 days continuous	6 pack	comfilcon A / 48%	Diameter	14.0mm
Sphere Power	Design		FDA Group	SiHy
+8.00D to -12.00D (0.50D steps after +/-6.00D)	Aspheric		UV Blocking*	No
+15.00D to -20.00D (0.50D steps after +/-6.00D)			Visibility Tint	Yes

Extraordinary visual performance. Premium all-day comfort.

Biofinity® toric



Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.7mm
Daily or up to 29 days continuous	3 pack	comfilcon A / 48%	Diameter	14.5mm
Sphere Power	Cylinder Power	Axis	FDA Group	SiHy
+8.00D to -10.00D (0.50D steps after +/-6.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	UV Blocking*	No
			Visibility Tint	Yes
			Design	Optimized Toric Lens Geometry™
+8.50 to +10.00 (0.50D steps)	-0.75 to -2.25 (in 0.50 steps)	5° to 180° (5° Steps)		
+10.00 to -10.00 (0.50D steps after +/- 6.00D)	-2.75 to -5.75 (in 0.50 steps)	5° to 180° (5° Steps)		

Extraordinary visual performance. Premium all-day comfort.

Biofinity® multifocal



Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.6mm
Daily wear or continuous wear up to 29 days	6 pack	comfilcon A / 48%	Diameter	14.0mm
Sphere Power	Add Power	Design	FDA Group	SiHy
+6.00D to -10.00D (0.50D steps after -6.00D)	+1.00 +1.50 +2.00 +2.50	D lens N lens	UV Blocking*	No
			Visibility Tint	Yes

Extraordinary visual performance at every distance. Premium all-day comfort.

Biofinity® toric multifocal




Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.6mm
Daily wear or continuous wear up to 29 days	3 pack	comfilcon A / 48%	Diameter	14.0mm
Sphere Power	Cyl Power	Design	FDA Group	SiHy
-10.00 to -6.50 (0.50 steps) -6.00 to +6.00 (0.25 steps) +6.50 to +10.00 (0.50 steps)	-0.75 to -5.75 (0.50 steps)	D lens N lens	UV Blocking*	No
	Add Power		Visibility Tint	Yes
	+1.00 +1.50 +2.00 +2.50			

Extraordinary visual performance. Premium all-day comfort.

Avaira Vitality™ family


Avaira Vitality™



High performance. Outstanding value.

Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.4mm
Monthly	6 pack	fanfilcon A / 55%	Diameter	14.2mm
Sphere Power	Design		FDA Group	SiHy
+08.00 to -12.00 (to +/- 6.00 in 0.25 steps) (from +/- 6.00 in 0.50 steps)	Aspheric front surface		UV Blocking*	Yes
			Visibility Tint	Yes

Avaira Vitality™ toric



High performance. Outstanding value.

Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.5mm
Monthly	3 pack	fanfilcon A / 55%	Diameter	14.5mm
Sphere Power	Cylinder Power	Axis	FDA Group	SiHy
-10.00D to +8.00D (0.50D steps after +/-6.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	UV Blocking*	Yes
		Design	Visibility Tint	Yes
		Optimized Toric Lens Geometry™		

Proclear® family

Proclear® 1 day



Clear vision*, convenience and comfort, for contact lens wearers who experience eye dryness.

Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.7mm
Daily	30 or 90 pack	omafilcon A / 60%	Diameter	14.2mm
Sphere Power	Design		FDA Group	2
+8.00D to -12.00D (0.50D steps after +5.00D and -6.00D)	Aspheric		UV Blocking*	No
			Visibility Tint	Yes

Proclear® 1 day multifocal



Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.7mm
Daily disposable	30 pack	omafilcon A / 60%	Diameter	14.2mm
Sphere Power	Add Power	Design	FDA Group	2
+6.00D to -10.00D (0.50D steps after -6.00D)	Single power profile, up to +2.50	Centre near aspheric	UV Blocking*	No
			Visibility Tint	Yes

Convenience and comfort for multifocal contact lens wearers who experience eye dryness.

Proclear®



Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.6mm
Daily	6 pack	omafilcon B / 62%	Diameter	14.2mm
Sphere Power	Design		FDA Group	2
+20.00D to -20.00D (0.50D steps after +/-6.00D)	Sphere		UV Blocking*	No
			Visibility Tint	Yes

Clear vision² and comfort for contact lens wearers who experience eye dryness.

Proclear® toric



Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.8mm
Daily	3 or 6 pack	omafilcon A / 62%	Diameter	14.4mm
Sphere Power	Cylinder Power	Axis	FDA Group	2
+6.00D to -8.00D (0.50D steps after -6.50D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	UV Blocking*	No
+6.25D to +10.00D (0.50D steps after +6.50D)	-0.75 to -2.25 (0.50D steps)	5° to 180° (5° steps)	Visibility Tint	Yes
-8.50D to -10.00D (0.50D steps)				
+10.00D to -10.00D (0.50D steps, +6.00D to -6.50D in 0.25D steps)	-2.75 to -5.75 (0.50D steps)	5° to 180° (5° steps)		

Helps provide all-day comfort for contact lens wearers with astigmatism who experience eye dryness.

Proclear® multifocal



Helps provide all-day comfort for multifocal contact lens wearers who experience eye dryness.

Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.7mm
Daily	3 or 6 pack	omafilcon A / 62%	Diameter	14.4mm
Sphere Power	Add Power	Design	FDA Group	2
+6.00 to -08.00	+1.00 +1.50 +2.00 +2.50	D lens N lens	UV Blocking*	No
+20.00D to -20.00D	+3.00 +3.50 +4.00		Visibility Tint	Yes
+6.00D to -6.00D (0.25D steps)				
+6.00D to +20.00D (0.50D steps)				
-6.50D to -20.00D (0.50D steps)				

Proclear® multifocal toric



Helps provide all-day comfort for multifocal contact lens wearers with astigmatism who experience eye dryness.

Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve	8.4, 8.8mm
Daily	3 pack	omafilcon B / 62%	Diameter	14.4mm
Sphere Power	Cylinder Power	Axis	FDA Group	2
+20.00D to -20.00D (0.50D steps after +/-6.50D)	-0.75 to -5.75 (0.50D steps)	5° to 180° (5° steps)	UV Blocking*	No
Design			Visibility Tint	Yes
D lens N lens				

* Due to aspheric optics.

¹ Brennan NA. Beyond Flux: Total Corneal Oxygen Consumption as an Index of Corneal Oxygen During Contact Lens Wear. Optom Vis Sci 2005. 2005;82(6):467-472.

² CVI data on file 2017 & 2014; non-dispensing short-term fitting study. N=40 & n=20.

* UV-blocking contact lenses help provide protection against transmission of harmful UV radiation to eye but are not substitutes for protective UV-absorbing eyewear, such as UV absorbing goggles or sunglasses, as they do not completely cover the eye or surrounding area. Continue to use UV-absorbing eyewear as directed by your eye care professional.